ARJOHUNTLEIGH

GETINGE GROUP





BOLERO™ - A RELIABLE LIFT BATH TROLLEY

Bolero makes assisted bathing routines for dependent residents a one-to-one care process.

For the resident, being bathed by a single carer is less stressful and intrusive. *Bolero* provides a secure environment throughout the bathing cycle, allowing the resident to relax and enjoy the full benefits of bathing. After collection from the bed, the resident remains safe on the lift bath trolley during transport, transfers and the bathing session itself.

A single carer can safely manage all stages of the bathing cycle. The efficient functions and features of *Bolero* ensure a safe, ergonomic working environment and give the carer more time to attend to the resident. With a paediatric version also available, *Bolero* can enhance the safety and comfort of bathing routines for dependent residents of all ages.

The resident is easily rolled onto the stretcher, which can be adjusted to the correct working height.



Bolero integrates with several height-adjustable ArjoHuntleigh bathing systems.



The stretcher can be fitted with a foldable side rail to enhance resident safety and comfort.



Handles on the pillar make the lift bath trolley easy to manoeuvre and position.



The quick-change battery system ensures Bolero is always available for use.



F

Electronic scale

Bolero can also be used to weigh residents using the optional foot-operated, highly accurate electronic scale

| Product Information | |
|----------------------|-----------------------------|
| Lifting capacity | 136 kg (300 lbs.) |
| Weight | 70 kg (154 lbs.) |
| Width | 825 mm (32 ½") |
| Length | 1850-1950 mm (72 ¾ - 76 ¾") |
| Length (paediatric) | 1350-1450 mm (53"-57") |
| Max total height | 1720 mm (67 ¾") |
| Lifting stroke range | 565 mm (22 ¼") |

| Powered functions | Raising & lowering |
|---------------------|---|
| Bath compatibility: | With Rhapsody/Primo (System 23), Rhapsody (System 25) and Classic Line |
| Best suited for | Mobility levels Doris and Emma |

Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.

@ and $^{\intercal M}$ are trademarks belonging to the ArjoHuntleigh group of companies. @ ArjoHuntleigh, 2016