
PRESSURE INJURY PREVENTION

ConformX®

Visco-elastic foam mattress range



A layer of visco-elastic foam provides support to the patient for pressure redistribution. The ConformX range was specifically developed with patient comfort and pressure injury management in mind, and provides a 24 hour solution.



ConformX Range

The ConformX range has been designed to provide a 24 hour care solution. The range includes a mattress replacement, mattress overlay and seat cushion. The ConformX range is available in various sizes and is easy to maintain due to the durability of the cover and foam.



Safe and convenient for patients and carers

Designed to provide effective pressure redistribution.

24 Hour Care Package Seat cushion available

The ConformX mattress replacement is available in multiple sizes and has been tested for X-Ray translucency¹ and is suitable for CPR². In addition, the overlay includes corner securing straps for extra security.

Cover features

Overlay System - A mattress overlay is

also available and is easy to install

- Ultra two way stretch
- Vapour permeable
- Water resistant
- Infection Control is aided by Antimicrobial properties in the cover material
- Non slip silver base
- Easy to clean and maintain
- Concealed zip
- Clearly marked covers with TurnTable[®] guide and washing instructions





TurnTable Guide - All ConformX mattress replacement systems carry the TurnTable Guide

2 CONFORMX 3

Specifications*

2-WAY TURN MATTRESS REPLACEMENT	DIMENSIONS H X W X D	ORDER CODE
MR Acute	1980 x 860 x 150 mm (78"x 33¾"x 6")	150CX198/086
MR Standard	1900 x 920 x 150 mm (74¾"x 36¼"x 6")	150CX190/092
MR Double	1900 x 1380 x 150 mm (74¾"x 54¼"x 6")	150CX190/138
MR Euro Wide	2040 x 920 x 150 mm (80¼"x 36¼"x 6")	150CX204/092
MR Euro Narrow	2000 x 800 x 150 mm (78¾"x 31½"x 6")	150CX200/080
Max. patient weight:	150 kg (330 lbs)	
Lifetime expectancy: 5 years for foam core and 2 years for cover		
Overlay systems		
OL Standard	1980 x 860 x 95 mm (78"x 33¾"x 3¾")	CO198/086
OL Community	1900 x 920 x 95 mm (78"x 33¾"x 3¾")	CO190/092
Seat cushions		
Seat	430 x 430 mm (17" x 17")	CS043/043
Max patient weight:	90 kg (198 lbs)	
Fire retardancy		

MR System meets BS EN 597 -1: 1995 BS EN 597 -2: 1995 and BS 6807: 1996, Section 2, Ignition Sources 0,1 and 5. Overlay System meets BS EN 597-1: 1995, BS EN 597 -2: 1995 and BS 6807: 1996, Section 2, Ignition Sources 0,1 and 5

Mattress

X-ray: Tested to FDA 21 CFR Ch. 1 (4-1-95 Edition) HHS & CEI/IEC 601-1-3: 1994-07 (X-ray attenuation). CPR: Chest Compression for CPR purposes is not impeded

References:

- $1. \ \, {\rm Arjo\,Data\,on\,File:}\, Conform X\,X ray\,Attenuation.\, Report\,Number\,100076011, April\,2002.$
- 2. Arjo Data on File Cardiopulmonary Resuscitation (CPR) for Foam Mattresses. Report Number 100076011, 2005

October 2019 Only Arjo designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by Arjo. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice. ® and M are trademarks belonging to the Arjo group of companies. © Arjo, 2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure injury and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

Arjo AB · Hans Michelsensgatan 10 · 211 20 Malmö · Sweden · +46 10 335 4500

www.arjo.com



 $^{{}^\}star\mathsf{Specifications}$ are subject to change without notice.