

ARJOHUNTLEIGH

GETINGE GROUP

MAXIONCE DISPOSABLE SLIDING AID




MaxiOnce™ is a single use disposable sliding aid designed to meet the needs and challenges of today's healthcare environments. Quality patient care is key – the same *MaxiOnce* is only used with one specific patient on one occasion, ensuring a hygienic and secure patient transfer. *MaxiOnce* is a tubular sliding aid made of a strong low-friction material. It is a great solution for lateral transfers in hospital areas with high patient turnover and high demand on infection control.

...with people in mind

www.ArjoHuntleigh.com

MAXIONCE – DISPOSABLE SLIDING AID



Accessibility & Availability

MaxiOnce comes on a roll, which is inserted into the *MaxiOnce* holder. The *MaxiOnce* can be cut to any length needed. Only the necessary length of material needs to be taken away from the roll before each transfer. The *MaxiOnce* measurement system makes it easy to determine the right length of the sliding aid.

For nursing and care staff, *MaxiOnce* is always available – no more wasted time looking for sliding aids. No more sliding aids will be damaged or lost during the laundry process. Most importantly, no more dangerous manual lifts because sliding aids are damaged or unavailable.

Reduces Cross-contamination

As single use disposables, *MaxiOnce* further protects against the risks associated with cross-contamination. Protecting patients is a key part of an infection control policy and further reduces the risk of patient litigation.

Cost Effective

Cost Management is high on the agenda for today's healthcare managers. *MaxiOnce* can help to provide safer patient transfers and reduce the risk of injuries to caregivers.



Between each ArjoHuntleigh logo, the distance is approximately 12 inches / 1 foot / 30 cm. The chart below shows the required length of *MaxiOnce* to be pulled from the roll.



The set-up enables an efficient use of the roll since only the necessary length of material needs to be taken from the roll before each transfer.



Article number: NMA1000



Article number: NMA2000

MaxiOnce	
Weight	8.4 kg (17.6 lbs.)
Dimension (Length x Width)	100 x 0.8 m (328 x 2½')

Length feet	Length meters	Number of ArjoHuntleigh logos
4 - 5	1.20 - 1.50	5
5 - 6	1.50 - 1.80	6
6 - 7	1.80 - 2.10	7
7 - 8	2.10 - 2.40	8

© and ™ are trademarks belonging to the ArjoHuntleigh group of companies.
© ArjoHuntleigh, 2009.